

PC Electrometer™

The Ultimate in Portable Reference Dosimetry

PC Electrometer is a dual channel reference class electrometer for absolute dose calibration. The system is designed for accuracy and convenience. It offers small size (0.4 kg), near no warm-up time (< 1 minute), and complete operation through USB, with no batteries or external power connections.



Features and Benefits


- Reference class dosimetry for absolute dose calibration
- Two independent measurement channels
- Lightweight and portable; only 0.4 kg
- USB powered – no batteries or power cord
- Intuitive user interface
- Interfaces with the 1D SCANNER™
- Less than 1-minute warm-up time
- Single USB cable connection
- Fast sampling interval of 500 ms
- Detector library

Automatic start/stop with threshold detection


Data display with configurable units and detector calibration library

Fully configurable interface and parameter options


Data logging capability




Features



PC Electrometer sampling is every 500 ms. This graph represents the dynamic range of the PC Electrometer when sampling is lowered to 100 ms (not available in commercial systems). Each step in this graph of total charge indicates a linac pulse as measured with a 0.6 cc chamber. PC Electrometer clearly shows the increase in linac pulse frequency as a function of dose rate. From 80 MU/min, the additional pulses are 1, 2, 3, and 4 for the respective rates of 160, 240, 320, 400.



Fast sampling makes PDD curves smooth and accurate when measured with the PC Electrometer and 1D SCANNER. The PDD curve above is an unsmoothed ratio of Field to Reference, using 0.6 cc Farmer type chambers.



1D SCANNER™



The PC Electrometer directly connects to the 1D SCANNER for monthly and annual QA.

1D SCANNER control

USB power/data when not using 1D SCANNER

8-Pin DIN connector. Used with PDI 2.0 or higher

Specifications

Warm Up Time: < 1.0 min

Charge Range: 2 pC – 10 mC, 15 fC resolution

Current Range (Continuous): Low: 2 pA – 50 nA

Current Range (Pulsed): 0 -105 pC/pulse

Leakage Drift: ± 0.001 pA

Display Update Frequency(s): 500 ms

Bias Voltage: Adjustable, 0 to ± 400 V

Non-linearity: $\pm 0.1\%$ of full scale

Long Term Stability: $\pm 0.5\%$

Measurement Repeatability: $\pm 0.25\%$ of full scale

A/D Converter: 16 bit

Operating System: Windows 10, Windows 8.1 or Windows 7 (32 or 64 bit)

Dimensions / Weight: 10.6 x 14.8 x 4.5 cm / 0.46 kg

Compatibility: SNC Dosimetry

Conformity: Reference class according to IEC 60731

PEO B.V.

info@gotopeo.com
www.gotopeo.com

The Netherlands

Havenweg 16, 6603 AS Wijchen
+31 (0)24 648 86 88

Belgium

Watermolenstraat 2, B-2910 Essen
+32 (0)3 309 32 09

CoC 34107894

VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46

BIC RABONL2U

