


SNC Machine™

Automate Your QA



AUTOMATE YOUR QA

SNC Machine listens for and captures your QA files, processes and analyzes the files, and saves the results to the database. Simply login to SNC Machine and immediately view a dashboard of results. Accept results that pass, and drill down into the analysis details for results that fail. Trend any piece of data against any other piece of data. It is that simple, and that powerful.



TG-142 & VMAT Imaging and Mechanical QA

All TG-142 recommended imaging and mechanical tests for monthly QA are included with SNC Machine. Simply deliver the test beam and SNC Machine does the rest. Accept or reject results on your terms and timeframe.

Tests Include

- Image Quality & Accuracy
 - CBCT
 - kV
 - MV
- MLC
 - Picket Fence
 - Log File Positioning
 - Leaf Speed
- Winston-Lutz Isocenter
 - Radiation
 - Machine
- Star shot
 - Gantry
 - Couch
 - Collimator
- Light / Radiation Field Congruence
- Beam
 - Field Size
 - Flatness
 - Symmetry

All VMAT tests recommended by Varian, but also applicable to Elekta, are included with SNC Machine. Simply deliver the test beam and SNC Machine does the rest. Accept or reject results on your terms and timeframe.

Tests Include

- Dose Rate versus Gantry Speed
- Leaf Speed
- Arc Point Dose
- DMLC Point Dose



19 different tests for TG-142 & VMAT QA




MLC Picket Fence





Star Shot

Results and images for a Star Shot Test with automatically calculated circle diameter and spoke angles.




Winston-Lutz

Results and images for a Winston-Lutz Isocenter Test. View the 2D results for each image as well as the combined 3D results of all images.




SNC Machine Trending

Any test results, such as CBCT, can be trended and visualized against other test results



CBCT

Results and images for a CBCT Test. Rollover any thumb nail to see the parameters and the specific corresponding ROI's to ensure image registration was performed correctly.



Trending

Trend any test parameter, against other test parameters, for any number of machines. View trends within the context of pass / fail criteria with base-lines and comments visible in the trend graph.

More SNC Machine Highlights

- Included SNC Server streamlines installation and provides access to any networked computer
- Easy phantom baseline setup with automatic ROI registration
- Supports most common imaging and mechanical QA phantoms for TG-142
- Works with Varian, Elekta, Aria, Mosaic

Specifications

Tests

TG-142 Imaging: CBCT Image Quality & Accuracy,
kV Image Quality & Accuracy,
MV Image Quality & Accuracy

TG-142 Mechanical: MLC Picket Fence, MLC Positioning & Leaf
Speed, Winston-Lutz Radiation & Machine
Isocenter, Gantry/Couch/Collimator
Starshot, Light/Radiation Congruence

Beam: Field Size, Beam Flatness, Beam Symmetry

VMAT: Dose Rate versus Gantry Speed, Leaf
Speed, Arc Point Dose, DMLC Point Dose

Phantoms

Sun Nuclear: MV-QA, kV-QA, FS-QA, WL-QA

Standard Imaging: PipsPro Phantoms


Phantom Laboratory: CatPhan 503, 504, 600

Leeds: TOR 18FG

Varian: Las Vegas Phantom

Gammex: 464

Sun Nuclear Phantoms




WL-QA

Dimensions: 60 x 60 x 60 mm

Sphere size: 8.0 mm

Sphere center accuracy: 20 mm





FS-QA

Field sizes: 100 x 100 mm;
150 x 150 mm

Markers (±0.1 mm): 56 - Field size (7 per field edge)
1 - Orientation

Dimensions: (L x W x D) 178 x 178 x 6 mm



MV-QA / kV-QA

MV Line pairs: 0.1, 0.2, 0.5, 1.0 ± 0.025 mm
kV Line pairs: 0.6, 1.2, 1.8, 2.4 ± 0.01 mm

MV ROI: 9 (4 spatial, 4 contrast, 1 center)
kV ROI: 28 (4 spatial, 23 contrast, 1 center)

MV Dimensions: (L x W x D) 127 x 102 x 25 mm

kV Dimensions: (L x W x D) 127 x 127 x 16 mm



radiation technology



SUN NUCLEAR
corporation

**YOUR MOST VALUABLE QA &
DOSIMETRY TOOLS**