

TomoDose™

TomoTherapy® Beam QA Made Easy


TomoDose is a 2D array for independent quality assurance and trouble-shooting measurements of a TomoTherapy® system. The system provides quick simultaneous full field profile measurements for the Y and X axes and comprehensive beam analysis.



Features and Benefits

- Use for QA on all TomoTherapy systems
 - Routine QA
 - Post component replacement QA
 - Daily QA
- Measure entire Hi-Art beam in a single measurement
 - Single X-axis and nine Y-axis detector arrays (total of 10)
 - Output, symmetry, flatness, field size, penumbra, beam center
 - Instantaneous dose rate
- Only array designed for TomoTherapy Hi-Art
- Accessibility to raw data
- Import water tank data and compare to TomoDose files
- Fast and efficient setup

TomoDose software displaying X and Y Axes



TomoDose_D052815

All data used is best available at time of publication. Data is subject to change without notice. ©2015 Sun Nuclear Corporation. All Rights Reserved.

PEO B.V.

info@gotopeo.com
www.gotopeo.com

The Netherlands

Havenweg 16, 6603 AS Wijchen
+31 (0)24 648 86 88

Belgium

Watmolenstraat 2, B-2910 Essen
+32 (0)3 309 32 09

CoC 34107894

VAT NL807859151B01


sunnuclear.com // +1 321 259 6862

Corporate Headquarters: 3275 Suntree Boulevard, Melbourne, FL 32940 USA

IBAN NL29 RABO 0356 1960 46
BIC RABONL2U



Features



Specifications

Detector Type:	SunPoint Diode Detectors
Detector Quantity:	223 total on X and Y
Detector Spacing: (mm)	X Primary: 5.0 Y Primary: 4.0 Y Offset: 8.0
Field Size:(cm)	53.0 x 9.8
Array Length: (cm)	X: 53.0 Y: 9.8 (8.0 at ± 19.0)
Y Axes Offset: (cm)	± 5.0 , ± 10.0 , ± 15.0 , ± 19.0 cm
Inherent Backscatter: (g/cm ²)	2.3
Active Detector Area: (mm ²)	0.64

Detector Volume: (cm ³)	0.000019
Detector Sensitivity: (nC/Gy)	32.0
Detector Stability:	0.5%/kGy at 6MV
Maximum Dose Rate Limit: (Gy/min)	56.0
Operating System:	Windows 2000, XP 32-bit, or Vista 32-bit
Dimensions: L/W/H (cm)	25.6 x 52.0 x 6.0
Weight: (kg)	5.0
Number of Connection Cables:	Single power/data cable

PEO B.V.

info@gotopeo.com
www.gotopeo.com

The Netherlands

Havenweg 16, 6603 AS Wijchen
+31 (0)24 648 86 88

Belgium

Watmolenstraat 2, B-2910 Essen
+32 (0)3 309 32 09

CoC 34107894
VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46
BIC RABONL2U

SunPoint®
PRECISION QA

