

GAMMEX

MAMMOGRAPHY



Artifact Identification Phantom

GAMMEX 179

The Gammex179 Artifact Identification Phantom is a full field device that allows a grey film to be produced without the phantom being imaged. This will indicate any artifacts that may be caused by an mammography system's grid or filters. It is also a great tool to verify that the film processor is working correctly. The Gammex179 is a great QA check tool to verify that your mammography system is working smoothly on all levels.

SPECIFICATIONS

Size..... 3.8x30x24cm (1.5x11.8x9.6 in)
 Weight..... 3.2kg (7 lbs)



WWW.GAMMEX.COM

PEO B.V.

info@gotopeo.com
 www.gotopeo.com

The Netherlands

Havenweg 16, 6603 AS Wijchen
 +31 (0)24 648 86 88

Belgium

Watermolenstraat 2, B-2910 Essen
 +32 (0)3 309 32 09

CoC 34107894

VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46

BIC RABONL2U

