

# Light field Ruler Set

Sun Nuclear (Gammex Technology)  
Beam geometry & alignment testing



Light Field Ruler Set measures the coincidence of the light and radiation fields of analog or digital flat plate/film X-ray units, including mammography systems.

Size/Range of Light Field to Measure	2 x 2 cm --> 25 x 30 cm with markings for 10, 15 x 15 & 20 x 20 cm and 25 x 25 cm
Size/Range of Circular Fields	1 cm – 25 cm diameter
Read Out Accuracy	0.5 mm
Dimensions	13.75 x 11 x 1/8 inches thick
Weight	938 g (2 lbs.)
Warranty	1 Year