



Isocentric Rotation Plate

GAMMEX 442-R

The Gammex 442-R Isocentric Rotation Plate is a quality assurance test instrument

designed to make necessary tests on radiotherapy machines quickly and easily. The Isocentric Rotation Plate takes only minutes to set up and performs routine QA tests such as:

- Light/Radiation Field Coincidence
- Symmetry and Flatness
- Accuracy of Optical Distance Indicators and Laser Alignment
- Gantry and Collimator Isocentricity

An optional ion chamber holder can be provided to be used with the base and side support for "in air" measurements.

SPECIFICATIONS

- Overall Size 45.7x31.8x34.9 cm
(18x12.5x13.75 in)
- Weight 4.4 kg (10 lbs)
- Ion Chamber
- Holder Optional RO5138 Inner diameter
1.27 cm (0.5 in)
- Film Used 10x12 in ready-pack