


GAMMEX



DIAGNOSTIC RADIOLOGY



# Processor QC Kit

## GAMMEX 185D

**Quality assurance in radiology begins with the film processor.** The processor is the single most influential source of problems in the diagnostic imaging department. To test all the parameters of the processor, Gammex provides the Gammex 185D Processor Quality Control Kit.

Included in this kit are a Portable Blue/Green Sensitometer, Portable Densitometer, Digital Thermometer, Fixer Retention Test kit and the text book "Film Processing in Medical Imaging". With these tools, daily processor quality control can be completed and early detection of possible problems can be caught before they lead to poor quality films or expensive repairs.

### SPECIFICATIONS

- The Gammex 185D includes:**
- Gammex MA0025 Densitometer
  - Gammex MA5034 Sensitometer
  - Gammex TM-99A Digital Thermometer
  - Gammex 081A Foam Lined Hard Carrying Case
  - Gammex 583 Book - "Film Processing in Medical Imaging"
  - Gammex 166B Fixer Retention Test Kit

Size . . . . . 44.5x18.4x35.6 cm  
(17.5x7.25x14 in)

Weight . . . . . Approx. 7 kg (15.5 lbs)

Optional Gammex 185D(HV) Kit contains a Portable Densitometer with 200 VAC charger.

