

GAMMEX



DIAGNOSTIC RADIOLOGY



Fluoroscopic Dose Rate and Low Contrast Resolution Test Tool Kit

GAMMEX 151

Comply with regulatory requirements for QA testing of fluoroscopic output and low contrast response with the Gammex 151 Low Contrast Resolution Test Tool Kit provides the test tools needed to comply with these requirements.

Low contrast performance in fluoroscopic units is crucial in making the right diagnosis in patients. Gammex has supplied the 151 Low contrast Resolution Test Tool to assure that your fluoroscopic system is working at maximum efficiency. Besides its elegant design, what sets the Gammex 151 apart from others is its ability to measure the automatic brightness stabilizer, maximum entrance exposure rate at the table top, and the phototimer performance of the spot film device.

Because of its ease of use and flexible design, our Low Contrast Resolution Test Tool can achieve accurate low contrast readings within the 2% range. By using the Gammex 151 in conjunction with our High Contrast Resolution Test Tool (Model 141 and 141H) you can insure that your fluoroscopic system is working at optimal efficiency and effectiveness. The components of the Gammex 151 Kit can also be used to test aspects of the performance of digital and analog spot-film devices, radiographic automatic exposure control functions, as well as digital angiography and digital cine-angiography.

SPECIFICATIONS

Construction Two Aluminum Blocks,
(1) Lead Blocker, (1) Aluminum
Resolution Plate
Dimensions 18x18x4.5 cm (7x7x1.8 in)
Weight 4 kg (8.8 lbs)

PEO B.V.

info@gotopeo.com
www.gotopeo.com

The Netherlands

Havenweg 16, 6603 AS Wijchen
+31 (0)24 648 86 88

Belgium

Watermolenstraat 2, B-2910 Essen
+32 (0)3 309 32 09

CoC 34107894
VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46
BIC RABONL2U

