

GAMMEX

Dual Energy CT Phantom Ca and I Contrast Quality Control

RADIATION ONCOLOGY



The Gammex Dual Energy CT Phantom provides users with the ability to perform Quality Assurance for Dual Energy CT analysis of Iodine and Calcium.

The Phantom consists of a Solid Water disk approximately the size of an average pelvis. A matrix of 16 holes in the disk hold interchangeable rods made of materials containing 7 different concentrations each of Iodine and Calcium. The rods can be positioned as the user chooses. The 2 remaining holes are occupied by Solid Water rods.

Scanning the phantom on a periodic basis provides data useful for the QA program related to the detect-

ability range of the Dual Energy CT scanner. The phantom uses the same base as the Gammex 467 Tissue Characterization Phantom, meaning owners of that phantom can take advantage of the interchangeability of the rods for the 2 products to provide a more economic approach to your Quality Assurance program

Gammex

WWW.GAMMEX.COM

PEO B.V.

info@gotopeo.com
www.gotopeo.com

The Netherlands

Havenweg 16, 6603 AS Wijchen
+31 (0)24 648 86 88

Belgium

Watermolenstraat 2, B-2910 Essen
+32 (0)3 309 32 09

CoC 34107894
VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46
BIC RABONL2U



GAMMEX 472

RADIATION ONCOLOGY



Calcium and Iodine rods of various concentrations permit easy analysis of detection ranges.

Specifications

Disk:

Material.....Gammex 451 Solid Water
Diameter.....33 cm (12.9 in)
Height.....5 cm (2 in)

Weight

Phantom + case.....6.6 kg (14.5 lbs)
Disk and Rods.....4.6 kg (10 lbs)

Rod Content

Iodine	Calcium
2.0 mg/ml	50.0 mg/ml
2.5 mg/ml	100.0 mg/ml
5.0 mg/ml	200.0 mg/ml
7.5 mg/ml	300.0 mg/ml
10.0 mg/ml	400.0 mg/ml
15.0 mg/ml	500.0 mg/ml
20.0 mg/ml	600.0 mg/ml

PEO B.V.

info@gotopeo.com
www.gotopeo.com

The Netherlands

Havenweg 16, 6603 AS Wijchen
+31 (0)24 648 86 88

Belgium

Watermolenstraat 2, B-2910 Essen
+32 (0)3 309 32 09

CoC 34107894
VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46
BIC RABONL2U

