

GAMMEX

CT Phantom Ring Extension

DIAGNOSTIC RADIOLOGY



Gammex 464-Ring - A simple inexpensive method of adding diagnostic function to your CT Quality Testing program

The Gammex 464 CT Phantom is a recognized standard for quality assurance monitoring of your CT scanner. Along with the optional analysis software, the phantom provides users with a wide range of quality benchmarks that can be measured.

The 464-Ring is an optional module that can be used with the phantom itself to provide a more anatomically correct image of scanning performance when focusing on the torso. This is achieved by extending the diameter of the traditional 464 phantom.

The Solid Water® composed Ring is designed to allow

the phantom to be easily inserted into it. All of the tests that can be performed on the phantom itself can also be performed when the optional Ring is used. The shape of the Ring is designed to mimic that of a normal human torso.

The Ring Extension, along with the 464 ACTS software, add further value to your 464 CT Phantom.



WWW.GAMMEX.COM

PEO B.V.

info@gotopeo.com
www.gotopeo.com

The Netherlands

Havenweg 16, 6603 AS Wijchen
+31 (0)24 648 86 88

Belgium

Watermolenstraat 2, B-2910 Essen
+32 (0)3 309 32 09

CoC 34107894

VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46

BIC RABONL2U



GAMMEX 464-Ring

DIAGNOSTIC RADIOLOGY



The 464 CT Phantom in the 464-Ring ready for scanning



The 464 CT Phantom being inserted into the 464-Ring

Extension Ring Dimensions and Specifics

Materials. Solid Water
 Dimensions
 Height. 10 3/8 in
 Width. 13 in
 Length. 6 3/8 in (cm)
 Weight. 13.5 lbs (kg)

Key Benefits of the Extension Ring

- Permits optional scanning with 464 phantom that better mimics a torso
- Easy to convert and use
- Compact and easy to store
- Economical