


Beta/ Gamma Laminate Credit Card Sources Laminated Sources

Spectrum Techniques

Laminate credit card sources offer a convenient alternative package which is easy to handle and store. Each source is constructed using 8 mil card stock and a laminate sleeve. The standard size is 3.75" x 2.25".

All beta/gamma laminate credit card sources have an uncertainty of ± 20% of the labeled activity. Activities will not exceed the U.S. NRC Exempt Quantity limit.

Not available for calibration



Available Nuclides and Activities

Barium 133
Cadmium 109
Cesium 137
Cobalt 57
Cobalt 60
Europium 152
Manganese 54
Sodium 22
Strontium 90
Thallium 204

Nuclide

Minimum Activity
0.1 µCi (3.7 kBq)
0.1 µCi (3.7 kBq)
0.05 µCi (1.85 kBq)
0.1 µCi (3.7 kBq)
0.05 µCi (1.85 kBq)
0.1 µCi (3.7 kBq)
0.1 µCi (3.7 kBq)
0.1 µCi (3.7 kBq)
0.05 µCi (1.85 kBq)
0.1 µCi (3.7 kBq)
0.1 µCi (3.7 kBq)

Maximum Activity
10 μCi (370 kBq)
10 μCi (370 kBq)
10 μCi (370 kBq)
100 μCi (3700 kBq)
1 µCi (37 kBq)
1 µCi (37 kBq)
10 μCi (370 kBq)
10 μCi (370 kBq)
0.1 µCi (3.7 kBq)
10 μCi (370 kBq)
10 µCi (370 kBq)

Zinc 65

+32 (0)3 309 32 09

