


# Mammography Phototimer Consistency Testing Slabs

Model 014A



## BETTER THAN PMMA FOR AEC CALIBRATION

The American Cancer Society and American College of Radiology guidelines for the screening of asymptomatic women have made over 50 million women candidates for mammography. In light of the staggering numbers involved, it's critically important that simple but reliable methods be developed to assess system performance and assure consistent production of diagnostically useful images.

CIRS Phototimer Consistency Testing Slabs are designed for precise assessment of AEC system performance in accordance with American College of Radiology and MQSA recommendations. Unlike acrylic, these testing slabs are manufactured with very tight thickness tolerances and more accurately simulate real breast tissue over the range of energies used in mammography. BR-12 (47% glandular / 53% adipose) is most commonly used.

### Features

- Assess AEC system performance
- Comply with ACR & MQSA recommendations
- Available in multiple configurations

900 Asbury Ave • Norfolk, Virginia 23513 • USA  
Tel: 800.617.1177 • 757.855.2765 • Fax: 757.857.0523

[WWW.CIRSINC.COM](http://WWW.CIRSINC.COM)

**CIRS**

Tissue Simulation & Phantom Technology

**PEO B.V.**

info@gotopeo.com  
www.gotopeo.com

**The Netherlands**

Havenweg 16, 6603 AS Wijchen  
+31 (0)24 648 86 88

**Belgium**

Watermolenstraat 2, B-2910 Essen  
+32 (0)3 309 32 09

CoC 34107894

VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46

BIC RABONL2U



## SPECIFICATIONS

### MAMMOGRAPHY PHOTOTIMER CONSISTENCY TESTING SLABS

MODEL NO.	STANDARD GLANDULARITY	QUANTITY	LENGTH (CM)	WIDTH (CM)	THICKNESS (CM)
014A	BR12	3	12.5	10	2
		2	12.5	10	1
		1	12.5	10	0.5

#### References:

White, D.R., R.J. Martin, and R. Darlison, Epoxy resin based tissue substitutes, British Journal of Radiology, 5, 814-821, 1977.

Materials are formulated to maximize simulation properties at 20 keV for the mammographic range, 80 keV for the diagnostic range and 0.5 MeV and above for the therapeutic range.

©2013 Computerized Imaging Reference Systems, Inc. All rights reserved.  
 Specifications subject to change without notice.  
 Publication: 014A DS120919



Computerized Imaging Reference Systems, Inc. has been certified by UL DQS Inc. to **(ISO) 13485:2016**. Certificate Registration No.10000905-MP2016..

#### PEO B.V.

info@gotopeo.com  
 www.gotopeo.com

#### The Netherlands

Havenweg 16, 6603 AS Wijchen  
 +31 (0)24 648 86 88

#### Belgium

Watermolenstraat 2, B-2910 Essen  
 +32 (0)3 309 32 09

CoC 34107894

VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46

BIC RABONL2U

