



Portable Blue/Green Sensitometer

GAMMEX MA5034

Sensitometry is the single most effective way to test the processor operation and consistency. For day to day processor quality control, the Gammex MA5034 Sensitometer provides a versatile, reliable, and highly accurate test. With its 21 step light modulator, a full range of densities can be tested with a single piece of film. The Gammex MA5034 is useful for all types of film, from sensitive x-ray film to roll or cine film. It can test either blue sensitive or green sensitive film. To expose the film simply close the cover and listen for the tone.

SPECIFICATIONS

Time Stability . . . ± 0.02 Log Exposure per Year
 Reproducibility . . . ± 0.04 Log Exposure
 Power 9 V Alkaline Battery (included)
 Warm-up Time . . None
 Blue Color Peak
 Wavelength 460 nm ± 10 nm
 Green Color Peak
 Wavelength 510 nm ± 10 nm
 Size 3.8x7.6x17.8 cm (1.5x3x7 in)
 Weight 0.57 kg (1.25 lbs)



WWW.GAMMEX.COM



GAMMEX MA5034

DIAGNOSTIC RADIOLOGY



0608 © Gammex, Inc. All rights reserved.

GAMMEX INC.
P.O. BOX 620327
MIDDLETON, WI 53562-0327
USA
+1 800 GAMMEX1 (426 6391)
+1 608 828 7000
FAX: +1 608 828 7500
EMAIL: SALES@GAMMEX.COM

GAMMEX-RMI LTD.
BROADWAY BUSINESS CENTRE
32A STONEY STREET
NOTTINGHAM NG1 1LL
UNITED KINGDOM
+44 (0) 115 924 7188
FAX: +44 (0) 115 924 7189
EMAIL: UKSALES@GAMMEX.COM

GAMMEX-RMI GMBH
FRANKFURTER STRASSE 15
D-35390 GIESSEN
GERMANY
+49 (0) 641 250 9176
FAX: +49 (0) 641 966 2642
EMAIL: DESALES@GAMMEX.COM

For more information contact GAMMEX today at 1 800 GAMMEX 1 (426 6391) or visit WWW.GAMMEX.COM

Nederland

T +31 (0)24 648 86 88

België / Belgique

T +32 (0)3 309 32 09

info@gotoPEO.com

www.gotoPEO.com