



Portable Densitometer

GAMMEX MA0025

Designed to meet today's demanding needs, the Gammex MA0025 Densitometer is portable, light weight, reliable and highly accurate. The Gammex MA0025 delivers impressive performance characteristics; a 0-4.0 optical density measuring range and an accuracy of ± 0.02 optical density. These rigid standards plus its rugged construction make the Gammex MA0025 ideal for field and hospital use alike. In addition, it comes with its own carrying case, certified stepwedge, and low battery indicator. A battery eliminator to accommodate AC power is also included.

SPECIFICATIONS

Range 0 to 4.0 Optical Density
 Accuracy ± 0.02 Density
 Reproducibility . . ± 0.01 Density
 Warm-up Time . . None
 Measuring Area . . 2 mm dia. and 1 mm dia.
 Power Supply . . . Four rechargeable AA Nicad Batteries, 4.8 V total rated at 600 mAh (included)
 Battery Charger . . SE 30-45 (115 VAC) or SE 30-46 (230 VAC) 50 to 60 Hz
 Charge Time Approx. 14 hrs.
 Size 5.08x7.46x17.8 cm (2x2.9x7 in)
 Weight 0.7 kg (1.5 lbs)

Certified Stepwedge included



WWW.GAMMEX.COM



GAMMEX MA0025

DIAGNOSTIC RADIOLOGY



0608 © Gammex, Inc. All rights reserved.

GAMMEX INC.
P.O. BOX 620327
MIDDLETON, WI 53562-0327
USA
+1 800 GAMMEX1 (426 6391)
+1 608 828 7000
FAX: +1 608 828 7500
EMAIL: SALES@GAMMEX.COM

GAMMEX-RMI LTD.
BROADWAY BUSINESS CENTRE
32A STONEY STREET
NOTTINGHAM NG1 1LL
UNITED KINGDOM
+44 (0) 115 924 7188
FAX: +44 (0) 115 924 7189
EMAIL: UKSALES@GAMMEX.COM

GAMMEX-RMI GMBH
FRANKFURTER STRASSE 15
D-35390 GIESSEN
GERMANY
+49 (0) 641 250 9176
FAX: +49 (0) 641 966 2642
EMAIL: DESALES@GAMMEX.COM

For more information contact GAMMEX today at 1 800 GAMMEX 1 (426 6391) or visit WWW.GAMMEX.COM

Nederland

T +31 (0)24 648 86 88

België / Belgique

T +32 (0)3 309 32 09

info@gotoPEO.com

www.gotoPEO.com